




Factory Address: 4# Workshop, No. 2875 Dongfu Avenue,  
Haicang District, Xiamen  
Tel.: +86 592 2293630  
Fax: +86 592 2296092  
eMail: sales@bellamoon.cn  
Website: <https://bellamoon.en.alibaba.com>

**倍乐梦（厦门）医疗科技有限公司**  
Bellamoon (Xiamen) Medical Technology Co.,Ltd



▶ Professional Manufacturer of Ward Nursing Equipment

倍乐梦（厦门）医疗科技有限公司  
Bellamoon (Xiamen) Medical Technology Co.,Ltd

## C1280J Medical gynecological examination bed

- Length of bed surface  $\geq 1420\text{mm}$
- Width of bed surface  $\geq 570\text{mm}$
- Lowest height  $\geq 600\text{mm}$
- Lift range of the bed surface  $\geq 600\text{mm}$
- Adjustment of leg-putting frame Random
- Power supply AC220V 50Hz
- Hand railing 1 Pair
- Leg-putting frame 1 Pair
- Dirt basin 1 Piece
- Waterproof gasket 1 Piece
- Power supply line 1 Piece
- Foot pedal 1 Piece
- Pedal switch 1 Piece



## C1280F Hospital multifunctional electric delivery bed

- Length of bed surface  $\geq 2100\text{mm}$
- Width of bed surface  $\geq 940\text{mm}$
- Lowest height  $\geq 660\text{mm}$
- Lift range of the bed surface  $\geq 300\text{mm}$
- Front inclination of the bed surface  $\geq 75^\circ$
- Back inclination of the bed surface  $\geq 5^\circ$
- Up-folding of a seat plate  $\geq 15^\circ$
- Up-folding of back plate  $\geq 15^\circ$
- Handle 1 Pair
- Leg-putting frame 1 Pair
- Dirt basin 1 Piece
- Rotating foot pedal 1 Pair
- Waterproof gasket 1 Piece
- Hand controller 1 Piece
- Power supply line 1 Piece
- Infusion support 1 Piece

## C1280Y Medical hydraulic gynecology chair

- Length of bed surface ≥1850mm
- Width of bed surface ≥600mm
- Length of auxiliary table ≥550mm
- Width of auxiliary table ≥500mm
- Lowest height of the bed surface ≥660mm
- Front inclination of the bed surface ≥10°
- Back inclination of the bed surface ≥15°
- Up folding of back board ≥70°
- Down folding of back board ≥18°
- Arm-putting frame 1 Pair
- Leg-putting frame 1 Pair
- Dirt basin 1 Piece
- Waterproof gasket 1 Piece



## C1280J2 Electric hydraulic gynecological operating table

- Length of bed surface ≥1850mm
- Width of bed surface ≥600mm
- Length of auxiliary table ≥550mm
- Width of auxiliary table ≥500mm
- Lowest height of the bed surface ≥600mm
- Lift range of the bed surface ≥330mm
- Auxiliary table expands outward ≥100°
- Seat plate adjustment -2°-10°
- Back board adjustment -6°-70°
- Arm-putting frame 1 Pair
- Leg-putting frame 1 Pair
- Handle 1 Pair
- Waterproof gasket 1 Piece
- Dirt basin 1 Piece
- Power supply line 1 Piece
- Pedal switch 1 Piece

## S1380F Medical portable surgical operating table

- Length of bed surface  $\geq 2100\text{mm}$
- Width of bed surface  $\geq 500\text{mm}$
- Lowest height of the bed surface  $\geq 650\text{mm}$
- Lift range of the bed surface  $\geq 250\text{mm}$
- Left inclination of bed surface  $\geq 20^\circ$
- Right inclination of the bed surface  $\geq 20^\circ$
- Front inclination of the bed surface  $\geq 12^\circ$
- Back inclination of the bed surface  $\geq 12^\circ$
- Reflex of head board  $\geq 50^\circ$
- Down folding leg plate  $\geq 90^\circ$
- Outward turning of leg plate  $\geq 90^\circ$
- Up folding of back board  $\geq 60^\circ$
- Lifting of waist bridge  $\geq 80\text{mm}$



## C1280Z6 Gynecological examination table with drawers

- Length of bed surface  $\geq 1850\text{mm}$
- Width of bed surface  $\geq 650\text{mm}$
- Bed surface height(customized)  $\geq 700\text{mm}$
- Up folding of back board  $\geq 60^\circ$
- Foot pedal 1 Pair
- Leg-putting frame 1 Pair
- Hand railing 1 Pair
- Auxiliary table 1 Piece
- Waterproof gasket 1 Piece