



Factory Address: 4# Workshop, No. 2875 Dongfu Avenue,
Haicang District, Xiamen
Tel.: +86 592 2293630
Fax: +86 592 2296092
eMail: sales@bellamoon.cn
Website: <https://bellamoon.en.alibaba.com>

倍乐梦（厦门）医疗科技有限公司
Bellamoon (Xiamen) Medical Technology Co.,Ltd

Product Catalogue V1.4

Bellamoon
倍乐梦



Professional Manufacturer of Ward Nursing Equipment

倍乐梦（厦门）医疗科技有限公司
Bellamoon (Xiamen) Medical Technology Co.,Ltd

TC5002A Medical ABS transfer treatment trolley

- Material:ABS Plastic, aluminum and steel
- The treatment cart uses 304 stainless steel guardrail, special concave top design, covered with transparent plastic.
- The car body is mainly composed of plastic, aluminum and steel structure, ABS double-layer flying swallow tabletop, novel design, beautiful appearance.
- Three drawers, one-piece ABS engineering material.
- Left side: a hand sanitizer basket is installed.
- Right side: Two sharp sharps boxes and two pedal trash cans are placed in the basket.
- The base is equipped with a large removable storage basket.



TC6001A Multifunctional ABS plastic medical cart

- Material: ABS plastic,Internal steel and aluminum structure.
- Guardrails on top of trolley is made of stainless steel, fix sharps box, avoid possible fall down during movement of trolley.
- With handrails on both side, convenient for pushing easily.
- Concave ABS moulding table, stainless steel guardrails, 2 sharp boxes, 1 side hanging box, 2 dirt barrels.
- 1 drawer, 1 middle workbench, one storage box below.
- Detachable storage box which is designed tightly fixed on the frame, and meanwhile easily take out whenever necessary.
- Foot-operated trash can with lid.



TC4001A Luxury ABS medical emergency trolley

- Material: Steel • ABS engineering plastic structure composition.
- The drawer is equipped with a one-time password lock.
- The cart is equipped with CPR board, suitable for emergency situations.
- Drawers: two small, two medium, and one large, each with compartments.
- Plug-in oxygen cylinder holder, strap fixed on the main body of the trolley.
- The infusion stand is inserted into the main column and can be interchanged with the position of the defibrillator.
- The emergency vehicle uses high-strength ABS material, and the plastic surface has strong resistance to corrosive chemicals.
- Medical emergency vehicles can store all types of first aid equipment and supplies, intravenous infusion stands, defibrillators, oxygen cylinders, trash cans, sharps boxes and medicines.



TC4001B Hospital ABS patient emergency cart

- Material: Aluminum • Steel • ABS engineering plastic structure composition;
- Size: 750*475*920mm(length * width*height)
- Drawer card holder: card insert from downside, easy keep clean.
- Plug-in oxygen cylinder holder, strap fixed on the main body of the trolley.
- The drawer is equipped with a one-time Code lock.
- The cart is equipped with CPR board, suitable for emergency situations.
- 5 drawers: two small, two medium, and one large drawer(inner with dividers can be organized freely into different size units)
- The infusion stand is inserted into the main column and can be interchanged with the position of the defibrillator.
- Aluminum columns with impact protection design and raised edges.



TC5001A Hospital medical nursing trolley

- Material: Aluminum • Steel • ABS engineering plastic structure composition.
- Size: 850*520*950mm (length * width * height)
- Nursing trolley use one-piece ABS plastic top board with raised-edge design, covered transparent soft plastic glass.
- IV pole, Fixed double-row 5-hook infusion rack.
- 2 small, 1 middle, 1 big drawers, inner with dividers can be organized easily and freely, label cards for the trolley and drawers.
- Equipped with four layers of drawers, built-in 3*3 separators, which can be separated freely, seal slot-type signage, to prevent liquid and dust entry.
- Centralized lock, stainless steel guard rail.
- Plastic basket, two dust baskets, hand washing liquid holder, Infusion rack, Sharps box.
- Four luxurious noiseless casters, two with brakes.



TC7002A Medical ABS hospital mobile medicine trolley

- Material: ABS body, aluminum structure, stainless steel guardrail.
- Size: 850*520*950mm (length * width * height)
- Medicine drawer: 25 small medicine boxes.
- Left side of car body: removable file box.
- Right side of the car body: double dirt buckets are used for garbage classification.
- The bottom of the car body: luxury universal plug-in mute wheels, two of which have brake function.
- Can be used in multiple departments: medicine delivery cart, treatment trolley.
- The upper part of the car body: ABS injection molding process table, with armrests on both sides, recessed design to prevent items from sliding off, stainless steel guardrail, and transparent soft glass on the table.



TC6005A Medical stainless steel equipment trolley

- Material:Electrostatic spraying on the surface of stainless steel.
- Size:690*450*940mm(length * width*height)
- Dressing trolley Color and size could be customized.
- 2 drawers with aluminum alloy handle.
- 2 plastic removable trash can.
- Stainless steel side rails prevent the trolley from falling when it moves. There is a stainless steel handrail on one side, which is convenient for carts.
- TPR tire no worn out after running 30KM, united forming without bolts,pass dynamic test: Bearing 120kg run 30KM and pass obstacles 500time.Casters are made of high-quality plastic material, light and flexible.



TC6009A Medical Stainless steel dressing trolley

- Material: ABS body, aluminum structure, stainless steel guardrail.
- Size:500*400*710/840mm(length * width*height)
- Dressing instrument trolley Color and size could be customized.
- 1 drawer with aluminum alloy handles.
- With PVC soft glass in drawer.
- Stainless steel side rails can prevent items from falling when the trolley is moving. Stainless steel guardrail is convenient for carts.
- TPR tire no worn out after running 30KM, united forming without bolts,pass dynamic test: Bearing 120kg run 30KM and pass obstacles 500time.Casters are made of high-quality plastic material, light and flexible.



TC9001A Hospital file patient records medical trolley

- Material: Body material made of ABS; Powder coated steel column bearing.
- Size: 750*410*980mm (length * width * height)
- Patient File Cart body is made of plastic and steel.
- With 2 drawers and 50 partitions for holding 50 clips, besides, with two lockers.
- Drawer easy pull out with high flexibility.
- According to user's demand to equip dossier.
- According to user's demand to increase or decrease the number of lattice.
- With key lock to prevent loss of case folder, used to lock all medical record folders.
- A striking number label side, convenient case management.



TC9001B Medical mobile plastic patient record trolley

- Material: ABS body, powder coated steel column bearing.
- Size: 450*410*950mm (length * width * height)
- With 1 drawers and 25 partitions for holding 25 folder, besides, with 1 lockers.
- Drawer easy pull out with high flexibility.
- According to user's demand to equip dossier.
- A striking number label side, convenient case management.
- Silent castors with brake.
- The top adopts ABS injection molding, with handrails to prevent items from turning over. There are handrails on both sides.



TC1105A waterproof canvas patient dirt cleaning trolley

- Material:Stainless steel
- Size:800*590*890mm(length * width*height)
- The hospital laundry trolley is made of powder coated steel material and designed without sharp edges.
- With two-piece waterproof canvas bag.
- A yellow bag and a blue bag, both equipped with lids.
- The top cover has strong tightness, anti-odor, ABS engineering material injection molding.
- Foot pedal opening and closing advanced cushioning design, closed silent.
- The cart is made of stainless steel, ABS engineering plastics, high-quality waterproof canvas, ultra-quiet casters.



TC1201A Hospital mobile cleaning linen trolley

- Material: Stainless steel.
- Size:588*428*890mm(length * width*height)
- Foot Lift Type Dirt Bag, Automatically Reset, Dirt Bag: Canvas, Anti-Leakage, Easy To Clean.
- Anti-Collision Devices Are Installed Around The Cart Body.
- Equipped With 4pcs Mute Castors and 2pcs Brake.
- Usage: Quilt & Clothes Cart, Garbage Cart, Syringe Cart, Dirt Cart.
- The Surface Is Easy To Clean and Wipe.



TC3001A Hospital multifunctional surgical anesthesia cart

- Material:aluminum, steel and ABS engineering plastic structure.
- Size:750*475*920mm(length * width*height)
- The anesthesia cart is composed of aluminum, steel and ABS engineering plastic structure; Aluminum alloy four-column load bearing.
- Front: Central lock, 5 drawers.
- Left side: Hidden retractable side table, detachable file box.
- Right side: a sharps box and 2 trash cans.
- Bottom: Luxurious medical silent wheels, two of which have brake function.
- ABS injection molding one-time molding, with handrails on both sides, stainless steel guardrail, transparent soft glass on the table.



TC6002A ABS steel multifunctional hospital treatment cart

- Material: ABS,stainless steel.
- Size:680*480*920mm(length * width*height)
- One-piece ABS plastic top board with raised-edge design, covered transparent soft plastic glass.
- 1 drawer with partitions inside, which can be easily organized.
- With four plastic-steel columns.
- 2*2 partitions inside the drawer, which can be separated freely.
- The tabbed area is designed according to the principles of ergonomics.



TC6003C Customized stainless steel trolley

- Material:Stainless steel.
- Size:690*450*940mm(length * width*height)
- Dressing trolley Color and size could be customized.
- 2 drawers with aluminum alloy handle.
- Plastic removable trash can.
- Stainless steel side rails prevent the trolley from falling when it moves. There is a stainless steel handrail on one side, which is convenient for carts.
- TPR tire no worn out after running 30KM, united forming without bolts,pass dynamic test: Bearing 120kg run 30KM and pass obstacles 500time.Casters are made of high-quality plastic material, light and flexible.



TC6004A stainless steel waterproof canvas nursing cart

- Material: 304 stainless steel
- Size:900*500*920mm(length * width*height)
- Made by stainless steel and fabric.
- Configure 4 φ75 medical silent casters, moves freely and easily.
- The specification could be made as your request.
- There are three stainless steel platform boards on the right and a waterproof canvas bag on the left, which is convenient and practical.
- The trolley is equipped with four truckles which can be braked respectively, and when in static state, it is stable.

